


Press Control



The complete monitoring and control tool
for hydrogenation type reactions

Features and Benefits

- Perform reactions at constant feed rate or constant pressure
- Use with any size/ style HP reactor
- Live gas uptake monitoring
- Simple single connection to your existing reactor



For more information or a demonstration please contact us
www.helgroup.com T: +44 (0) 20 8441 6778 E: marketing@helgroup.com

Press Control Includes


- Pressure control unit with mass flow sensor
- Temperature and pressure sensor
- Electronics and PC with monitoring and control software

Ease of use


- Connect your gas feeds to the control unit
- Connect the single outlet line to your reactor
- Perform your reaction while recording all the data

Data handling

- Real time data display of gas uptake
- The data is automatically saved to a data file
- The data is readily exported to Excel for analysis



**Press Control Unit with PC monitoring
and control software**



Typical data trace from Press Control hydrogenation

For more information or a demonstration please contact us
www.helgroup.com T: +44 (0) 20 8441 6778 E: marketing@helgroup.com