




6050 Compensated Jacket Calorimeter



The new Model 6050 Compensated Jacket Calorimeter is the first calorimeter designed by Parr to be controlled by your own personal computer, helping to make it Parr's most affordable calorimeter to date.

Features of the new 6050 Compensated Jacket Calorimeter include:

- PC Control: The 6050 is attached directly to a Windows PC using a standard USB cable. All functions and tests are accessible with the click of a mouse using Parr's 6050 Software Package, installed on your own PC
- Automatic Oxygen Filling: The 6050 fills the combustion vessel with the correct amount of oxygen for each test automatically – removing one manual step from the standard combustion vessel preparation routine
- Operates at room temperature
- Operator time per test: 6 minutes

Let Us Build One for You.